

TAIPEI MEDICAL UNIVERSITY

CORE FACILITY CENTER, OFFICE OF RESEARCH AND DEVELOPMENT

Inductively Coupled Plasma Mass Spectrometer

AUTHENTICATION RECORD

APPLICANT PROFILE

NAME:		DATE:	
DEPARTMENT:		TEL. NUMBER:	
EMAIL:		ADVISOR SIGNATURE:	

AUTHENTICATION PROCEDURE:

1. Complete the training courses.
 - A. Join live, record the date and time. Then signed by CFC Staff.
 - B. Watch on-demand video on website (TMU account will be needed for sign in). Be sure to record the date and time, then signed by your advisor.
2. Pass the written examination.
3. Accumulate 1–5 operating experiences.
4. Pass the operating test.
5. Present an authentication record and apply for authorization.

PLEASE FOLLOW THESE RULES:

1. To qualify for authentication: complete the training courses.
2. Please provide authentication within 1 year of completing the training courses.
3. Accumulate 1–5 operating experiences supervised by technicians or authorized operators of the instrument.
4. Charges apply during the accumulation of operational experience and the operating test.
5. Please accept supervision until passing the operating test.
6. Operating the instrument without authorization and supervision shall result in a 1 month loss of operating authorization for all involved parties.

Authentication items	Date and time	Sig. of supervisor
Training courses		
Written examination		
Operating experiences		
Operating test		

Approbation of Core Facility Center _____